



Transcatheter closure of VSD (my experiences)

**Nguyen Lan Hieu, MD, PhD
Hanoi Medical University
Binh Duong General Hospital**

Key points

- The transcatheter occlusion of VSD is considered to be one of the most sophisticated and complex interventional procedures.
- Complications of this procedure are Aortic regurgitation and AVB
- The small to moderate size of VSD specially in Adult patient is very common in Asia country

DIFFERENT DEVICES: DIFFERENT RESULTS?



PERIMEMBRANOUS VSD

3 MAIN POINTS



- Membranous septal aneurysm
- Distance to aortic valve
- Size of defect

PRESENCE OF ANEURYSM

Without aneurysm:


- Single opening exit in to RV
- ADO I, ADO II, Pfm coil, symmetric or asymmetric VSDO
- **Easy to cross the defect by guide wire or catheter**

WITHOUT MSA + SINGLE OPENING EXIT



WITHOUT MSA + SINGLE OPENING EXIT

ADO I, ADO II, PFM coil




Im: 2/77
Se: 0001
Lossy compression (JPEG)



WL: 128 WW: 256 [D]

PHAM HUU THIEN DUY
PHAM 24_24494
21-Oct-12
BV NHI DONG I
123454





21-Aug-14 12:34:54 PM

6

1

	
2 L 0.6 cm	
1 L 1.0 cm	










PRESENCE OF ANEURYSM

- Some space for the device
- Separating from aortic valve
- Tricuspid valve involvement in to MSA
- Triangle shape: PDA- like: ADO I, ADO II, coil, symmetric, asymmetric...
- Simple or complicated opening exits


WITH ANEURYSM




Some space for the device



ANEURYSM










FOR COMPLICATED ANEURYSM:

CLOSE ALL IN THE LEFT SIDE

MULTIPLE OPENING EXITS







MULTIPLE OPENING EXITS




TRICUSPID VALVE

- THE LENGTH OF ADO I : 7 – 8 mm
- BECOMES LONGER WHEN COMPRESSED








SMALL DEFECT AND ANEURYSM

- Difficult for guidewire crossing
- Difficult for big sheath
- Small and flexible profile is suitable
- Suitable for ADO II

Im: 2/77
Se: 0001
Lossy compression (JPEG)

DANG THUY NGAN
DANG 17_17653
25-Nov-03
BV NHI DONG I
104053




-80

WL: 128 WW: 256 [D]

21-Aug-14 10:40:53 AM

SMALL DEFECT



Distance from the aortic valve


❖ > 3 mm : many devices can be used: ADO I,
PFM COIL, ADO II, AMVSDO, Symmetric

Im: 2/77
Se: 0001
Lossy compression (JPEG)



WL: 128 WW: 256 [D]


PHAM HUU THIEN DUY
PHAM 24_24494
21-Oct-12
BV NHI DONG I
123454



~60

21-Aug-14 12:34:54 PM

ENOUGH AORTIC RIM




Distance from the aortic valve

❖ < 3 mm : need the soft and flexible device for
avoiding aortic injury: ADO II and COIL,
AAPMVSDO

Im: 2/77
Se: 0001
Lossy compression (JPEG)



NGO LE TU QUYNH
NGO L16_16312
12-Mar-14 0F
BV NHI DONG I
101832




WL: 128 WW: 256 [D]

21-Aug-14 10:18:32 AM



SUB-AORTIC: DEFICIENT AORTIC RIM



SUB-AORTIC: NO AORTIC RIM







NO AORTIC RIM




SUB-AORTIC: NO AORTIC RIM







RISK OF AORTIC REGURGITATION



Post MI VSD



VENTRICULAR SEPTAL RUPTURE COMPLICATING POST AMI





FOV 20
LAO 83


Residual VSD Post Operation



Post -Op: LV RA communication









Future direction for VSD closure

- **Sub Aortic or Pulmonary VSD**
- **New device with low profile**
- **International registry for VSD closure**

Softer and new design VSD device

Cocoon Device for Sub pulmonary VSD





Early results of new device

Results	Pts (n = 28)	%
Successfully implanted	27	96,42
Failure	1	2,58
Severe complication	1 (ischemic stroke Total recovery)	
Mild complication	1 PAC	

Follow up

	Before procedure (n=28)	After procedure (n=27)	6 months (n=24)	12 months (n=11)
Shunt		48,15 %	12,5%	9,1%
Minor AR	5 (17,86%)	4 (14,81%)	5 (20,83%)	2 (18,2%)
Moderate AR	0	1 (3,7%)	0	0
Arrhythmia	No	No	No	No

CONCLUSIONS

- **NO IDEAL DEVICE FOR VSD CLOSURE**
- **WIDE VARIETY OF VSD MORPHOLOGY**
- **MAIN KEY:**

DEVICE SELECTION AND PATIENT SELECTION

THANK YOU